Menicon
BLOOM[™]

Livre blanc

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* Bloom, Bloom Day, Bloom Night, Bloom Care, Menicon Z and Easyfit are trademarks of Menicon Co., Ltd.

Références

1. Holden BA, Fricke TR, Wilson DA, et al. Global prevalence of myopia and high myopia and temporal trends from 2000 through 2050. *Ophthalmology* 2016;123(5):1036-42.
2. Flitcroft DI. The complex interactions of retinal, optical and environmental factors in myopia aetiology. *Prog Retin Eye Res.* 2012;31(6):622-60.
3. Tano Y. Pathologic myopia: where are we now? *Am J Ophthalmol* 2002;134(5):645-60.
4. Vongphanit J, Mitchell P, Wang JJ. Prevalence and progression of myopic retinopathy in an older population. *Ophthalmology* 2002;109(4):704-11.
5. Wong TY, Klein BE, Klein R, Knudtson M, Lee KE. Refractive errors, intraocular pressure, and glaucoma in a white population. *Ophthalmology* 2003;110(1):211-7.
6. Saw SM, Gazzard G, Shih-Yen EC, Chua WH. Myopia and associated pathological complications. *Ophthalmic Physiol Opt* 2005;25(5):381-91.
7. Ikuno Y, Jo Y, Hamasaki T, Tano Y. Ocular risk factors for choroidal neovascularization in pathologic myopia. *Invest Ophthalmol Vis Sci* 2010;51(7):3721-5.
8. Tideman JW, Snabel MCC, Tedja MS, et al. Association of axial length with risk of uncorrectable visual impairment for europeans with myopia. *JAMA Ophthalmol* 2016; 134(12):1355-63.
9. Flitcroft DI, He M, Jonas JB, et al. IMI - Defining and classifying myopia: a proposed set of standards for clinical and epidemiologic studies. *Invest Ophthalmol Vis Sci* 2019;60(3):M20-M30.
10. Gifford KL, Richdale K, Kang P, et al. IMI - Clinical Management Guidelines Report. *Invest Ophthalmol Vis Sci* 2019;60(3):M184-M203.
11. Tedja MS, Haarman AEG, Meester-Smoor MA, et al. IMI - Myopia genetics report. *Invest Ophthalmol Vis Sci* 2019;60(3):M89-M105.
12. Troilo D, Smith 3rd EL, Nickla DL, et al. IMI - Report on experimental models of emmetropization and myopia. *Invest Ophthalmol Vis Sci* 2019;60(3):M31-M88.
13. Wolffsohn JS, Kollbaum PS, Berntsen DA, et al. IMI – Clinical myopia control trials and instrumentation report. *Investig Ophthalmol Vis Sci* 2019; 60(3):M132-M160.
14. Jones L, Drobe B, González-Méjome JM, et al. IMI – Industry guidelines and ethical considerations for Myopia control report. *Investig Ophthalmol Vis Sci* 2019;60(3):M161-M183.
15. Wildsoet CF, Chia A, Cho P, et al. IMI - Interventions for controlling myopia onset and progression report. *Invest Ophthalmol Vis Sci* 2019;60(3):M106-M131.
16. Goss DA, Winkler RL. Progression of myopia in youth: age of cessation. *Am J Optom Physiol Opt* 1983;60(8):651-8.
17. Group C. Myopia stabilization and associated factors among participants in the Correction of Myopia Evaluation Trial (COMET). *Invest Ophthalmol Vis Sci* 2013;54(13):7871-84.
18. Flaxman SR, Bourne RRA, Resnikoff S, et al. Global causes of blindness and distance vision impairment 1990-2020: a systematic review and meta-analysis. *Lancet Glob Heal* 2017;5(12):e1221-e34.
19. Vitale S, Cotch MF, Sperduto R, Ellwein L. Costs of refractive correction of distance vision impairment in the United States, 1999-2002. *Ophthalmology* 2006;113(12):2163-70.
20. Lim MC, Gazzard G, Sim EL, Tong L, Saw SM. Direct costs of myopia in Singapore. *Eye* 2009;23(5):1086-9.
21. Pan CW, Dirani M, Cheng CY, Wong TY, Saw SM. The age-specific prevalence of myopia in Asia: a meta-analysis. *Optom Vis Sci* 2015;92(3):258-66.
22. Koh V, Yang A, Saw SM, et al. Differences in prevalence of refractive errors in young Asian males in Singapore between 1996-1997 and 2009-2010. *Ophthalmic Epidemiol* 2014;21(4):247-55.
23. Pan CW, Ramamurthy D, Saw SM. Worldwide prevalence and risk factors for myopia. *Ophthalmic Physiol Opt* 2012;32(1):3-16.
24. Wang TJ, Chiang TH, Wang TH, Lin LL, Shih YF. Changes of the ocular refraction among freshmen in National Taiwan University between 1988 and 2005. *Eye* 2009;23(5):1168-9.
25. Vitale S, Sperduto RD, Ferris 3rd FL. Increased prevalence of myopia in the United States between 1971-1972 and 1999-2004. *Arch Ophthalmol* 2009;127(12):1632-9.

26. Gilmartin B. Myopia: precedents for research in the twenty-first century. *Clin Exp Ophthalmol* 2004;32(3):305-24.
27. Rose KA, Morgan IG, Smith W, Burlutsky G, Mitchell P, Saw SM. Myopia, lifestyle, and schooling in students of Chinese ethnicity in Singapore and Sydney. *Arch Ophthalmol* 2008;126(4):527-30.
28. French AN, Morgan IG, Mitchell P, Rose KA. Risk factors for incident myopia in Australian schoolchildren: the Sydney adolescent vascular and eye study. *Ophthalmology* 2013;120(10):2100-8.
29. Mutti DO, Hayes JR, Mitchell GL, et al. Refractive error, axial length, and relative peripheral refractive error before and after the onset of myopia. *Invest Ophthalmol Vis Sci* 2007;48(6):2510-9.
30. He M, Zheng Y, Xiang F. Prevalence of myopia in urban and rural children in mainland China. *Optom Vis Sci* 2009;86(1):40-4.
31. Gwiazda J, Hyman L, Dong LM, et al. Factors associated with high myopia after 7 years of follow-up in the Correction of Myopia Evaluation Trial (COMET) Cohort. *Ophthalmic Epidemiol* 2007;14(4):230-7.
32. Pacella R, McLellan J, Grice K, Del Bono EA, Wiggs JL, Gwiazda JE. Role of genetic factors in the etiology of juvenile-onset myopia based on a longitudinal study of refractive error. *Optom Vis Sci* 1999;76(6):381-6.
33. Wu MM, Edwards MH. The effect of having myopic parents: an analysis of myopia in three generations. *Optom Vis Sci* 1999;76(6):387-92.
34. Rosenfield M, Gilmartin B. *Myopia and Nearwork*. Oxford, UK.: Butterworth-Heinemann; 1998.
35. Rose K, Harper R, Tromans C, et al. Quality of life in myopia. *Br J Ophthalmol* 2000;84(9):1031-4.
36. Huang J, Wen D, Wang Q, et al. Efficacy comparison of 16 interventions for myopia control in children: a network meta-analysis. *Ophthalmology* 2016;123(4):697-708.
37. Smith 3rd EL. Optical treatment strategies to slow myopia progression: effects of the visual extent of the optical treatment zone. *Exp Eye Res* 2013;114:77-88.
38. Wen D, Huang J, Chen H, et al. Efficacy and acceptability of orthokeratology for slowing myopic progression in children: a systematic review and meta-analysis. *J Ophthalmol* 2015;2015:360806.
39. Li SM, Kang MT, Wu SS, et al. Studies using concentric ring bifocal and peripheral add multifocal contact lenses to slow myopia progression in school-aged children: a meta-analysis. *Ophthalmic Physiol Opt* 2017;37(1):51-9.
40. Hair LA, Steffensen EM, Benoit JS, Berntsen DA. Changes in peripheral defocus with center-distance and center-near multifocal contact lenses. *Optom Vis Sci* 2019:E-abstract 195209.
41. Lee R, Achenbach P, Schwiegerling J. Instantaneous power profiles vs sagittal power profiles of various center distance multifocal lenses. *Optom Vis Sci* 2019;96:E-abstract 195249.
42. Nti AN, Gregory HR, Ritchley ER, Berntsen DA. Contrast sensitivity with center-distance multifocal contact lenses. *Optom Vis Sci* 2019;96:E-abstract 190033.
43. Lopes-Ferreira D, Ribeiro C, Neves H, et al. Peripheral refraction with dominant design multifocal contact lenses in young myopes. *J Optom* 2013;6(2):85-94
44. Ladage PM, Yamamoto K, Ren DH, et al. Effects of rigid and soft contact lens daily wear on corneal epithelium, tear lactate dehydrogenase, and bacterial binding to exfoliated epithelial cells. *Ophthalmology* 2001;108(7):1279-88.
45. Morgan PB, Maldonado-Codina C, Efron N. Comfort response to rigid and soft hyper-transmissible contact lenses used for continuous wear. *Eye Contact Lens* 2003;29(1 Suppl):S127-30;
46. Maldonado-Codina C, Morgan PB, Efron N, Efron S. Comparative clinical performance of rigid versus soft hyper Dk contact lenses used for continuous wear. *Optom Vis Sci* 2005;82(6):536-48.
47. Morgan PB, Efron N, Maldonado-Codina C, Efron S. Adverse events and discontinuations with rigid and soft hyper Dk contact lenses used for continuous wear. *Optom Vis Sci* 2005;82(6):528-35.
48. Albright RA, Venuti BD, Ichijima H, Nyunt AK, Cavanagh HD. Postmarket surveillance of Menicon Z rigid gas-permeable contact lenses for up to 30 days continuous wear in the United States. *Eye Contact Lens* 2010;36(5):241-4.
49. Swarbrick HA. Orthokeratology review and update. *Clin Exp Optom* 2006;89(3):124-43.

50. Boost M, Cho P, Lai S. Efficacy of multipurpose solutions for rigid gas permeable lenses. *Ophthalmic Physiol Opt* 2006;26(5):468-75.
51. Choy CKM, Cho P, Boost M V. Cytotoxicity of rigid gas-permeable lens care solutions. *Clin Exp Optom* 2013;96(5):467-71.
52. Shi G Sen, Boost M V, Cho P. Does the presence of QAC genes in staphylococci affect the efficacy of disinfecting solutions used by orthokeratology lens wearers? *Br J Ophthalmol* 2016;100(5):708-12.
53. Cooper J, O'Connor B, Watanabe R, et al. Case Series Analysis of Myopic Progression Control With a Unique Extended Depth of Focus Multifocal Contact Lens. *Eye Contact Lens* 2018;44(5):e16-e24.
54. R. Payor, J. Woods, D. Fonn, P. Situ, S. Dillehay, R. Griffin, M. Tyson LJ. Feasibility Testing of a Novel Soft Contact Lens Optical Design to Reduce Suspected Risk Factors for the Progression of Juvenile Onset Myopia. *Invest Ophthalmol Vis Sci* 2014;55:E-abstract 3638.
55. S. Dillehay, J. Woods, P. Situ, S. Guthrie, R. Paynor, R. Griffin, M. Tyson, L. Jones. Comparison of Three Power Levels of A Novel Soft Contact Lens Optical Design to Reduce Suspected Risk Factors for the Progression of Juvenile Onset Myopia. *Invest Ophthalmol Vis Sci* 2014;55:E-abstract 3637.
56. Miller J, Long B, Dillehay S. Children's evaluation of a unique myopia progression control lens design. *Optom Vis Sci* 2013;88:E-abstract 115896.
57. Cho P, Cheung SW. Retardation of myopia in Orthokeratology (ROMIO) study: a 2-year randomized clinical trial. *Invest Ophthalmol Vis Sci* 2012;53(11):7077-85.
58. Santodomingo-Rubido J, Villa-Collar C, Gilmartin B, Gutiérrez-Ortega R. Myopia control with orthokeratology contact lenses in Spain: Refractive and biometric changes. *Invest Ophthalmol Vis Sci* 2012;53(8):5060-5
59. Chen C, Cheung SW, Cho P. Myopia control using toric orthokeratology (TO-SEE study). *Invest Ophthalmol Vis Sci* 2013;54(10):6510-7.
60. Santodomingo-Rubido J, Villa-Collar C, Gilmartin B, Gutiérrez-Ortega R, Sugimoto K. Long-term efficacy of orthokeratology contact lens wear in controlling the progression of childhood myopia. *Curr Eye Res* 2017;42(5):713-20
61. Bullimore MA, Sinnott LT, Jones-Jordan LA. The risk of microbial keratitis with overnight corneal reshaping lenses. *Optom Vis Sci* 2013;90(9):937-44.
62. Li SM, Kang MT, Wu SS, et al. Efficacy, safety and acceptability of orthokeratology on slowing axial elongation in myopic children by meta-analysis. *Curr Eye Res* 2016;41(5):600-8.
63. Liu YM, Xie P. The Safety of Orthokeratology--A Systematic Review. *Eye Contact Lens* 2016;42(1):35-42.
64. Bullimore MA. The Safety of Soft Contact Lenses in Children. *Optom Vis Sci* 2017;94(6):638-46.
65. Brennan NA, Toubouti Y, Greenaway N, Cheng X. Safety of contact lenses in pediatrics: data analysis from six randomized myopia control trials. *Optom Vis Sci* 2019;96:E-abstract: 195208.
66. Santodomingo-Rubido J, Villa-Collar C, Gilmartin B, Gutiérrez-Ortega R. Orthokeratology vs. spectacles: adverse events and discontinuations. *Optom Vis Sci* 2012;89:1133-9
67. Clinical Evaluation Report for Menicon Z Night Orthokeratology Contact Lenses. 2017.
68. Perrigin J, Perrigin D, Quintero S, Grosvenor T. Silicone-acrylate contact lenses for myopia control: 3-year results. *Optom Vis Sci* 1990;67(10):764-9.
69. Katz J, Schein OD, Levy B, et al. A randomized trial of rigid gas permeable contact lenses to reduce progression of children's myopia. *Am J Ophthalmol* 2003;136(1):82-90.
70. Walline JJ, Gaume A, Jones LA, et al. Benefits of contact lens wear for children and teens. *Eye Contact Lens* 2007;33(6):317-21.
71. Lipson MJ. Long-term clinical outcomes for overnight corneal reshaping in children and adults. *Eye Contact Lens* 2008;34(2):94-9.
72. Walline JJ, Jones LA, Sinnott L, et al. Randomized trial of the effect of contact lens wear on self-perception in children. *Optom Vis Sci* 2009;86(3):222-32.
73. Jones-Jordan LA, Chitkara M, Coffey B, et al. A comparison of spectacle and contact lens wearing times in the ACHIEVE study. *Clin Exp Optom* 2010;93(3):157-63.

74. Rah MJ, Walline JJ, Jones-Jordan LA, et al. Vision specific quality of life of pediatric contact lens wearers. *Optom Vis Sci* 2010;87(8):560-6.
75. Walline JJ, Jones LA, Sinnott LT. Corneal reshaping and myopia progression. *Br J Ophthalmol* 2009;93(9):1181-5.
76. Hiraoka T, Kakita T, Okamoto F, Takahashi H, Oshika T. Long-term effect of overnight orthokeratology on axial length elongation in childhood myopia: a 5-year follow-up study. *Invest Ophthalmol Vis Sci* 2012;53(7):3913-9.
77. Downie LE, Lowe R. Corneal reshaping influences myopic prescription stability (CRIMPS): an analysis of the effect of orthokeratology on childhood myopic refractive stability. *Eye Contact Lens* 2013;39(4):303-10.
78. Hiraoka T, Kakita T, Okamoto F, Oshika T. Influence of ocular wavefront aberrations on axial length elongation in myopic children treated with overnight orthokeratology. *Ophthalmology* 2015;122(1):93-100.
79. Lee YC, Wang JH, Chiu CJ. Effect of Orthokeratology on myopia progression: twelve-year results of a retrospective cohort study. *BMC Ophthalmol* 2017;17(1):243.
80. Hiraoka T, Sekine Y, Okamoto F, Mihashi T, Oshika T. Safety and efficacy following 10-years of overnight orthokeratology for myopia control. *Ophthalmic Physiol Opt* 2018;38(3):281-9.
81. Efron N, Morgan PB, Woods CA, International contact lens prescribing survey consortium. Survey of contact lens prescribing to infants, children, and teenagers. *Optom Vis Sci* 2011;88(4):461-8.
82. Wolffsohn JS, Calossi A, Cho P, et al. Global trends in myopia management attitudes and strategies in clinical practice. *Contact Lens Anterior Eye* 2016;39(2):106-16.
83. Efron N, Morgan PB, Woods CA, Santodomingo-Rubido J, Nichols JJ. International survey of contact lens fitting for myopia control in children. *Contact Lens Anterior Eye* 2019 [Epub ahead of print]
84. Wolffsohn J, Calossi A, Cho P, et al. Global trends in myopia management attitudes and strategies in clinical practice - 2019 Update. *Contact Lens Anterior Eye* 2019 [Epub ahead of print]
85. Santodomingo-Rubido J, Villa-Collar C, Gilmartin B, Gutiérrez-Ortega R. Myopia control with orthokeratology contact lenses in Spain: A comparison of vision-related quality-of-life measures between orthokeratology contact lenses and single-vision spectacles. *Eye Contact Lens* 2013;39(2):153-7
86. Bullimore M, Ritchey ER. Myopia control: an evidence-based comparison of the benefits and the risks. *Optom Vis Sci*. 2019;96:E-abstract 190031.
87. Mutti DO, Sinnott LT, Mitchell GL, et al. Relative peripheral refractive error and the risk of onset and progression of myopia in children. *Invest Ophthalmol Vis Sci*.2011;52(1):199-205.
88. Zadnik K, Sinnott LT, Cotter SA, et al. Prediction of Juvenile-Onset Myopia. *JAMA Ophthalmol*. 2015;133(6):683-689.
89. Donovan L, Sankaridurg P, Ho A, Naduvilath T, Smith 3rd EL, Holden BA. Myopia progression rates in urban children wearing single-vision spectacles. *Optom Vis Sci* 2012;89(1):27-32.
90. Cho P, Cheung SW. Discontinuation of orthokeratology on eyeball elongation (DOEE). *Contact Lens Anterior Eye* 2017;40(2):82-7.
91. Company Data on File, 2015.
92. Chan KY, Cheung SW, Cho P. Clinical performance of an orthokeratology lens fitted with the aid of a computer software in Chinese children. *Cont Lens Anterior Eye* 2012;35(4):180-4.
93. Wittich CM, Burkle CM, Lanier WL. Ten common questions (and their answers) about off-label drug use. *Mayo Clin Proc* 2012;87(10):982-90.
94. Aronson JK, Ferner RE. Unlicensed and off-label uses of medicines: definitions and clarification of terminology. *Br J Clin Pharmacol* 2017; 83(12):2615-25.
95. Riley Jr. JB, Basilius PA. Physicians' liability for off-label prescriptions. *Nephrol News Issues* 2007;21(7):43-4,46-7.
96. Wilkes M, Johns M. Informed consent and shared decision-making: a requirement to disclose to patients off-label prescriptions. *PLoS Med* 2008;5(11):e223.
97. Lenk C, Koch P, Zappel H, Wiesemann C. Off-label, off-limits? Parental awareness and attitudes towards off-label use in paediatrics. *Eur J Pediatr* 2009;168(12):1473-8.

98. Lenk C, Duttge G. Ethical and legal framework and regulation for off-label use: European perspective. *Ther Clin Risk Manag* 2014;10:537-46.
99. Weda Hoebert J, Vervloet M, Moltó Puigmarti C, Damen N, Marchange S, Langedijk J, Lisman J, van Dijk LM. Study on off-Label use of medicinal products in the European Union. © European Union 2017. Available at: https://ec.europa.eu/health/sites/health/files/files/documents/2017_02_28_final_study_report_on_off-label_use_.pdf. Accessed: December 19, 2019.

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S.A.S AU CAPITAL DE 12 523 900 EUROS - RCS BOBIGNY 722 038 668 - CODE INTRACOMMUNAUTAIRE FR92722038668